



Golden Vac™ Vacuum Blood Collection System

As a leading manufacturer of medical and laboratory plastic consumables in China over 36 years with ISO13485:2016 and CE certification, Gongdong continuously invests in new and innovated technologies to improve quality control and customers satisfaction throughout whole procedures.

- Fully automated production lines distribute to reliable and accurate digital inspection individually using a robotic imaging system
- Automatic labeling machines provide speed and high-precision by intelligent control with optimal readability Provide customized services for different requirements as various dimensions, labeling, or precise adjustments compatible with increasing altitude
- End-user and distributor training plans provided by our professional team
- Centrifugation:
≤ 1300 g (RCF) for 10 minutes at 18-25° C.



VACUUM BLOOD COLLECTION SYSTEM



Golden Vac™ Closures



- Various stopper colors conform with international standards.
- Rubber stopper is recessed inside the plastic shield to avoid potential contact from any blood drops left by a collection needle.
- High quality rubber stopper provides maximum sealing with easy puncturing and reliable self-sealing after needle piercing and preventing any possible interactions with blood samples in tubes.
- With lower residuals and impurities by suitable rubber formulation, rubber stopper chips seldom occur during needle piercing. And it also provides lower insertion forces and lower pulling forces.
- Gongdong patented non screw safety twist caps with strong, durable construction can be removed with a gentle pull, making outer plastic shield and interior rubber stopper more secured.
- Universal closures/safe-lock closures/rubber stoppers are available.
- Compatible with multi-types of automated analyzers.



Golden Vac™ Tubes



- Versatile tubes and reliable additives for optimal performance in the laboratory.
- Tubes are made of medical grade virgin polyethylene terephthalate (PET), clear as glass with equal or better performance, offer low permeability to solvents, eliminating breakage when tubes are centrifuged at higher speed or when stopper is removed and reinserted, to reduce risk of tube breakage during handling.
- Interior tube wall is coated with microscopic silica particles which activate clotting when tubes are gently inverted.
- Suitable for incineration or autoclaving; no harmful gas is generated.



Golden Vac™ Separator Gel

- High quality and reliable inert separation gel is ideal for the collection, storage, and transport of undiluted plasma for molecular diagnostics testing. When the tube is centrifuged, the separation gel migrates and forms a stable barrier between the plasma and most of the cellular elements without any interference to blood testing.

Preset Vacuum & High-Quality Additives

- Preset standard ration additives based on CLSI and CCCLS recommendation.
- Accurate preset vacuum level in the tubes ensure precise blood drawing.
- Various additives are ultrasonically spray-dried on interior wall to make sure coagulation or anticoagulation more thoroughly.



Vacuum Tubes for High Altitude

- Customized vacuum tube for high altitude to make sure accurate blood drawing in different area.





Serum Clot Activator Tube

Golden Vac™ serum clot activator tube is to obtain serum specimen used for serum determinations in chemistry, serology, and immunochemistry.

Description

- Silica particles act as the clot activator.
- Blood clotting promoted with ultrasonically spray-dried silica particles on the interior tube wall.
- Gently mix the blood collection tube by inverting 5-6 times immediately to activate clotting after blood drawing.
- To obtain serum, tube should be allowed to stand for 30 minutes or until the clot is formed.
- **Centrifugation:** ≤ 1300 g (RCF) for 10 minutes at 18-25° C.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|---------|----------------|-----------|-----------|-----------|-------------------|
| GD020CA | 2.0 | 13X75 | Clot Activator | / | PET/GLASS | | Safety Cap+Rubber |
| GD030CA | 3.0 | 13X75 | Clot Activator | / | PET/GLASS | | Safety Cap+Rubber |
| GD040CA | 4.0 | 13X75 | Clot Activator | / | PET/GLASS | | Safety Cap+Rubber |
| GD050CA | 5.0 | 13X75 | Clot Activator | / | PET/GLASS | | Safety Cap+Rubber |
| GD060CA | 6.0 | 13X100 | Clot Activator | / | PET/GLASS | | Safety Cap+Rubber |
| GD070CA | 7.0 | 13X100 | Clot Activator | / | PET/GLASS | | Safety Cap+Rubber |
| GD090CA | 9.0 | 16X100 | Clot Activator | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100CA | 10.0 | 16X100 | Clot Activator | / | PET/GLASS | | Safety Cap+Rubber |
| GD090CAR | 9.0 | 16X100 | Clot Activator | / | PET/GLASS | | Rubber Stopper |
| GA0100CAR | 10.0 | 16X100 | Clot Activator | / | PET/GLASS | | Rubber Stopper |

- Available in a range of 1-10ml blood collection tubes.
- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

Plain Tube

Golden Vac™ plain tubes with no additives for biochemistry, immunology and serology test.

Description

- Naturally separate out the blood cells and fibrin.
- The intact serum is suitable for a wide range of lab operation.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|----------|-----------|---------|----------|-----------|-----------|-----------|-------------------|
| GD020A | 2.0 | 13X75 | / | / | PET/GLASS | | Safety Cap+Rubber |
| GD030A | 3.0 | 13X75 | / | / | PET/GLASS | | Safety Cap+Rubber |
| GD040A | 4.0 | 13X75 | / | / | PET/GLASS | | Safety Cap+Rubber |
| GD050A | 5.0 | 13X75 | / | / | PET/GLASS | | Safety Cap+Rubber |
| GD060A | 6.0 | 13X100 | / | / | PET/GLASS | | Safety Cap+Rubber |
| GD070A | 7.0 | 13X100 | / | / | PET/GLASS | | Safety Cap+Rubber |
| GD090A | 9.0 | 16X100 | / | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100A | 10.0 | 16X100 | / | / | PET/GLASS | | Safety Cap+Rubber |
| GD090AR | 9.0 | 16X100 | / | / | PET/GLASS | | Rubber Stopper |
| GA0100AR | 10.0 | 16X100 | / | / | PET/GLASS | | Rubber Stopper |

- Available in a range of 1-10ml blood collection tubes.
- All tubes are supplied in 1000pcs or 1200 pcs/ctn.



EDTA Tube

Golden Vac™ EDTA Tubes contain EDTA as an anticoagulant, are widely used for determinations in EDTA whole blood clinical hematology test, especially for routine blood examination.

Description

- The EDTA inhibits the coagulation process by eliminating the calcium in the blood.
- Available in K2 · EDTA or K3 · EDTA.
- The anticoagulant effect of K3 · EDTA is stronger than that of K2 · EDTA. K3 · EDTA results in greater RBC shrinkage and produces a larger increase in cell volume on standing.
- Ultrasonically and finely coated to make sure accurate additives.
- After the tube has been filled with blood, immediately invert the tube 8-10 times to mix and ensure adequate anticoagulation of the specimen.
- **Centrifugation:** 1100 g for 10 minutes at 18-25° C (Only for EDTA tubes with separator gel)



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|------------|-----------|---------|-----------|-----------|-----------|-----------|-------------------|
| GD020EK2 | 2.0 | 13X75 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD030EK2 | 3.0 | 13X75 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD040EK2 | 4.0 | 13X75 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD050EK2 | 5.0 | 13X75 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD060EK2 | 6.0 | 13X100 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD070EK2 | 7.0 | 13X100 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD080EK2 | 8.0 | 16X100 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD090EK2 | 9.0 | 16X100 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100EK2 | 10.0 | 16X100 | K2 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100EK2R | 10.0 | 16X100 | K2 · EDTA | / | PET/GLASS | | Rubber Stopper |

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|------------|-----------|---------|-----------|-----------|-----------|-----------|-------------------|
| GD020EK3 | 2.0 | 13X75 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD030EK3 | 3.0 | 13X75 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD040EK3 | 4.0 | 13X75 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD050EK3 | 5.0 | 13X75 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD060EK3 | 6.0 | 13X100 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD070EK3 | 7.0 | 13X100 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD080EK3 | 8.0 | 16X100 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD090EK3 | 9.0 | 16X100 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100EK3 | 10.0 | 16X100 | K3 · EDTA | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100EK3R | 10.0 | 16X100 | K3 · EDTA | / | PET/GLASS | | Rubber Stopper |

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|------------|-----------|---------|-----------|-----------|-----------|-----------|-------------------|
| GD035EK2G | 3.5 | 13X75 | K2 · EDTA | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD060EK2G | 6.0 | 13X100 | K2 · EDTA | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD085EK2G | 8.5 | 16X100 | K2 · EDTA | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD085EK2RG | 8.5 | 16X100 | K2 · EDTA | Gel | PET/GLASS | | Rubber Stopper |
| GD035EK3G | 3.5 | 13X75 | K3 · EDTA | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD060EK3G | 6.0 | 13X100 | K3 · EDTA | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD085EK3G | 8.5 | 16X100 | K3 · EDTA | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD085EK3RG | 8.5 | 16X100 | K3 · EDTA | Gel | PET/GLASS | | Rubber Stopper |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

VACUUM BLOOD COLLECTION SYSTEM




Serum Separation Tubes

Golden Vac™ Gold-top vacuum blood tubes for serum determinations in chemistry, immunology and serology.

Description

- Coated with ultrasonically spray-dried microscopic silica particles which activate blood clotting.
- Contain an inert polymer gel that separates the serum from the bloodclot during centrifugation, preventing contamination of the serum.
- Gently mix the blood collection tube by inverting 5-6 times immediately to activate clotting after blood drawing.
- Allow the tube to stand for 30 minutes or until the clot is formed before centrifuging.
- **Centrifugation:** 1800-2200 g for 8-15 minutes at 18-25° C.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|---------|----------------|-----------|-----------|-----------|-------------------|
| GD040SGC | 4.0 | 13X75 | Clot Activator | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD040SGC2 | 4.0 | 13X100 | Clot Activator | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD050SGC | 5.0 | 13X100 | Clot Activator | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD060SGC | 6.0 | 13X100 | Clot Activator | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD080SGC | 8.0 | 16X100 | Clot Activator | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD080SGCR | 8.0 | 16X100 | Clot Activator | Gel | PET/GLASS | | Rubber Stopper |
| GD085SGC | 8.5 | 16X100 | Clot Activator | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD085SGCR | 8.5 | 16X100 | Clot Activator | Gel | PET/GLASS | | Rubber Stopper |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

ESR Tube

Determination of the Erythrocyte Sedimentation Rate for fully automated analyzers. ESR is a clinical lab test that measures the rate at which red blood cells in whole blood descend into a standardized tube, reported as mm per hour, to measure inflammation as inflamed red blood cells settle quickly relative to normal blood cells. A high ESR could be indicative of elevated inflammation.

Description

- Sodium citrate is used as anticoagulant to prevent blood from clotting.
- Filled with a buffered sodium citrate solution in a concentration of 0.129mol/L (3.8%).
- The blood to additive ratio is 4:1.
- After ESR tube has been filled with blood, immediately invert the tube 8-10 times to mix and ensure adequate anticoagulation of the specimen.
- Aluminum composite bags available to provide optimum protection against additive evaporation, UV light, dirt and mechanical influences.
- All the necessary custom service options to fulfill your project requirements.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|---------|-----------|---------|---------------------|-----------|-----------|-----------|-------------------|
| GD016SC | 1.6 | 13X75 | Sodium Citrate 3.8% | / | PET/GLASS | | Safety Cap+Rubber |
| GD024SC | 2.4 | 13X75 | Sodium Citrate 3.8% | / | PET/GLASS | | Safety Cap+Rubber |
| GD032SC | 3.2 | 13X75 | Sodium Citrate 3.8% | / | PET/GLASS | | Safety Cap+Rubber |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

VACUUM BLOOD COLLECTION SYSTEM



ESR Tube

For the manual determination of erythrocyte sedimentation rates with ESR stand.

Description

- Sodium citrate is used as anticoagulant to prevent blood from clotting.
- Filled with a buffered sodium citrate solution in a concentration of 0.129mol/L (3.8%).
- The blood to additive ratio is 4:1.
- After ESR tube has been filled with blood, immediately invert the tube 8-10 times to mix and ensure adequate anticoagulation of the specimen.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|---------|---------------------|-----------|----------|-----------|-------------------|
| GD016ESR | 1.6 | 8X120 | Sodium Citrate 3.8% | / | GLASS | | Safety Cap+Rubber |
| GD0128ESR | 1.28 | 8X120 | Sodium Citrate 3.8% | / | GLASS | | Safety Cap+Rubber |
| GD016ESRN | 1.6 | 9X120 | Sodium Citrate 3.8% | / | GLASS | | Safety Cap+Rubber |


- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

ESR Tube Rack

Immediately before the tubes are placed in the ESR tube rack for measurement, the tubes must be mixed again. After one or two hours the results are read. The ESR tube rack holds a maximum of 10 ESR tubes and has a height adjustable zero mark. The measurement results achieved correspond to the Westergren method.

| Cat. No | Size/mm |
|---------|---|
| H1020A | ESR Tube Rack, for ESR 8X120mm, 1.6ml tube |
| H1020B | ESR Tube Rack, for ESR 8X120mm, 1.28ml tube |

- All tube racks are supplied in cases of 10pcs.



Coagulation Tube

Golden Vac™ blue-top blood collection tubes are used for clinical coagulation tests.

Description

- Sodium citrate is used as anticoagulant to prevent blood from clotting.
- Filled with a buffered sodium citrate solution in a concentration of 0.109mol/L (3.2%).
- The blood to additive ratio is 9:1. The ratio of blood to anticoagulant is critical for valid prothrombin time results.
- Immediately after draw, gently invert the tube 3-4 times in order to activate adequate anticoagulant.
- Aluminum composite bags available to provide optimum protection against additive evaporation, UV light, dirt and mechanical influences.
- All the necessary custom service options to fulfill your project requirements.
- **Centrifugation:** 2000-2500g for 10-15 minutes at 18-25° C.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|----------|---------------------|-----------|-----------|-----------|-------------------|
| GD018SC | 1.8 | 13X75 | Sodium Citrate 3.2% | / | PET/GLASS | | Safety Cap+Rubber |
| GD027SC | 2.7 | 13X75 | Sodium Citrate 3.2% | / | PET/GLASS | | Safety Cap+Rubber |
| GD036SC | 3.6 | 13X75 | Sodium Citrate 3.2% | / | PET/GLASS | | Safety Cap+Rubber |
| GD045SC | 4.5 | 13X75 | Sodium Citrate 3.2% | / | PET/GLASS | | Safety Cap+Rubber |
| GD018SC-F | 1.8 | 10.25X64 | Sodium Citrate 3.2% | / | PET/GLASS | | Safety Cap+Rubber |
| GD027SC-F | 2.7 | 10.25X64 | Sodium Citrate 3.2% | / | PET/GLASS | | Safety Cap+Rubber |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.
- All tubes (Flat Bottom) are supplied in 1000pcs/ctn.



Double-wall tubes

Golden Vac™ blue-top blood tubes specially designed with double walls are made of PP (polypropylene) and PET (polyethylene terephthalate). PP is ideal for sensitive coagulation parameters due to its nonreactive characteristics. The inner tube is made of PP and prevents the evaporation of liquid. The outer tube is made of PET and ensures a long shelf-life for the vacuum.

Description

- Sodium citrate is used as anticoagulant to prevent blood from clotting. Filled with a buffered sodium citrate solution in a concentration of 0.109mol/L (3.2%).
- The blood to additive ratio is 9:1. The ratio of blood to anticoagulant is critical for valid prothrombin time results.
- Immediately after draw, gently invert the tube 3-4 times in order to activate adequate anticoagulant.
- **Centrifugation:** 2000-2500 g for 10-15 minutes at 18-25° C.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|----------|-----------|---------|---------------------|-----------|-------------------|-----------|-------------------|
| GD018SCD | 1.8 | 13X75 | Sodium Citrate 3.2% | / | PET (Double Wall) | | Safety Cap+Rubber |
| GD027SCD | 2.7 | 13X75 | Sodium Citrate 3.2% | / | PET (Double Wall) | | Safety Cap+Rubber |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

Heparin Tube

Golden Vac™ heparin tube is used for drawing heparinized plasma or whole blood for emergency biochemical and hemorheological tests.

Description

- Available in lithium heparin and sodium heparin.
- The anticoagulant heparin activates antithrombin, which blocks the clotting cascade and thus produces a whole blood/plasma sample.
- Lithium heparin has no interference with all ion detection during testing, including sodium ion.
- Immediately after draw, gently invert the tube 8-10 times in order to activate adequate anticoagulant.
- **Centrifugation:** ≤ 1300 g (RCF) for 10 minutes at 18-25° C.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|---------|-----------------|-----------|-----------|-----------|-------------------|
| GD020LH | 2.0 | 13X75 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD030LH | 3.0 | 13X75 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD040LH | 4.0 | 13X75 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD050LH | 5.0 | 13X75 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD060LH | 6.0 | 13X100 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD070LH | 7.0 | 13X100 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD080LH | 8.0 | 16X100 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD090LH | 9.0 | 16X100 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100LH | 10.0 | 16X100 | Lithium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100LHR | 10.0 | 16X100 | Lithium Heparin | / | PET/GLASS | | Rubber Stopper |

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|---------|-----------|---------|----------------|-----------|-----------|-----------|-------------------|
| GD020SH | 2.0 | 13X75 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD030SH | 3.0 | 13X75 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD040SH | 4.0 | 13X75 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD050SH | 5.0 | 13X75 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD060SH | 6.0 | 13X100 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD070SH | 7.0 | 13X100 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD080SH | 8.0 | 16X100 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |

VACUUM BLOOD COLLECTION SYSTEM



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|---------|----------------|-----------|-----------|-----------|-------------------|
| GD090SH | 9.0 | 16X100 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100SH | 10.0 | 16X100 | Sodium Heparin | / | PET/GLASS | | Safety Cap+Rubber |
| GD0100SHR | 10.0 | 16X100 | Sodium Heparin | / | PET/GLASS | | Rubber Stopper |

- Available in a range of 1-10ml blood collection tubes
- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

Heparin Tube (Separator Gel)

These Golden Vac™ heparin tubes contain a polymer separation gel and spray-dried heparin additive.

Description

- Available in lithium heparin and sodium heparin.
- The anticoagulant heparin activates antithrombin, which blocks the clotting cascade and thus produces a whole blood/plasma sample.
- Lithium heparin has no interference with all ion detection during testing, including sodium ion.
- The separation gel separates the plasma from the blood cells during centrifugation, preventing contamination of the plasma.
- Immediately after draw, gently invert the tube 8-10 times in order to activate adequate anticoagulant.
- **Centrifugation:** 1800-2200g for 8-15 minutes at 18-25° C.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|----------|-----------|---------|-----------------|-----------|-----------|-----------|-------------------|
| GD030LHG | 3.0 | 13X75 | Heparin Lithium | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD040LHG | 4.0 | 13X75 | Heparin Lithium | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD050LHG | 5.0 | 13X100 | Heparin Lithium | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD085LHG | 8.5 | 16X100 | Heparin Lithium | Gel | PET/GLASS | | Safety Cap+Rubber |

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|----------|-----------|---------|----------------|-----------|-----------|-----------|-------------------|
| GD030SHG | 3.0 | 13X75 | Sodium Heparin | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD040SHG | 4.0 | 13X75 | Sodium Heparin | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD050SHG | 5.0 | 13X100 | Sodium Heparin | Gel | PET/GLASS | | Safety Cap+Rubber |
| GD085SHG | 8.5 | 16X100 | Sodium Heparin | Gel | PET/GLASS | | Safety Cap+Rubber |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

Glucose Tube

Golden Vac™ glucose tube is used for blood glucose, glucose tolerance, red blood cell electrophoresis, sugar hemoglobin, etc.

Description

- Obtain highly stable plasma for blood glucose testing.
- Sodium Fluoride is an inhibitor of glucose degradation.
- Potassium Oxalate, EDTA Na₂, and Heparin Sodium prevent blood coagulation and avoid blood hemolysis.
- Immediately after draw, gently invert the tube 8-10 times in order to activate adequate anticoagulant.
- **Centrifugation:** ≤ 1300 g (RCF) for 10 minutes at 18-25° C.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|----------|-----------|---------|----------------------------------|-----------|-----------|-----------|-------------------|
| GD020SFH | 2.0 | 13X75 | Sodium Fluoride + Heparin Sodium | / | PET/GLASS | | Safety Cap+Rubber |
| GD030SFH | 3.0 | 13X75 | Sodium Fluoride + Heparin Sodium | / | PET/GLASS | | Safety Cap+Rubber |
| GD040SFH | 4.0 | 13X75 | Sodium Fluoride + Heparin Sodium | / | PET/GLASS | | Safety Cap+Rubber |
| GD050SFH | 5.0 | 13X75 | Sodium Fluoride + Heparin Sodium | / | PET/GLASS | | Safety Cap+Rubber |

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|----------|-----------|---------|--|-----------|-----------|-----------|-------------------|
| GD020SNA | 2.0 | 13X75 | Sodium Fluoride + EDTA Na ₂ | / | PET/GLASS | | Safety Cap+Rubber |
| GD030SNA | 3.0 | 13X75 | Sodium Fluoride + EDTA Na ₂ | / | PET/GLASS | | Safety Cap+Rubber |
| GD050SNA | 5.0 | 13X75 | Sodium Fluoride + EDTA Na ₂ | / | PET/GLASS | | Safety Cap+Rubber |

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|---------|-----------|---------|-------------------------------------|-----------|-----------|-----------|-------------------|
| GD020SP | 2.0 | 13X75 | Sodium Fluoride + Potassium Oxalate | / | PET/GLASS | | Safety Cap+Rubber |
| GD030SP | 3.0 | 13X75 | Sodium Fluoride + Potassium Oxalate | / | PET/GLASS | | Safety Cap+Rubber |
| GD040SP | 4.0 | 13X75 | Sodium Fluoride + Potassium Oxalate | / | PET/GLASS | | Safety Cap+Rubber |
| GD050SP | 5.0 | 13X75 | Sodium Fluoride + Potassium Oxalate | / | PET/GLASS | | Safety Cap+Rubber |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

VACUUM BLOOD COLLECTION SYSTEM



Non Vacuum Blood Tube

Micro collection tube and Capillary blood collection tube belongs to Non-Vacuum blood tube.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|---------|------------------|-----------|----------|-----------|-----------|
| GDN050EK2 | 5.0 | 18X55 | EDTA K2 | / | PP | | Screw Cap |
| GDN050EK3 | 5.0 | 18X55 | EDTA K3 | / | PP | | Screw Cap |
| GDN050A | 5.0 | 18X55 | Plain | / | PP | | Screw Cap |
| GDN050CA | 5.0 | 18X55 | Clot Activator | / | PP | | Screw Cap |
| GDN050SP | 5.0 | 18X55 | Fluoride Oxalate | / | PP | | Screw Cap |
| GDN050LH | 5.0 | 18X55 | Heparin Lithium | / | PP | | Screw Cap |

Non-Vacuum Serum Separator Tube

For serum determinations in chemistry, immunology and serology.

Description

- Contain an inert polymer gel that separates the serum from the blood clot during centrifugation, preventing contamination of the serum.
- Made of transparent polypropylene for easy visibility of blood samples PET.
- Immediately after draw, gently invert the tube 8-10 times in order to activate adequate anticoagulant.
- To obtain serum, tube should be allowed to stand for 30 minutes or until the clot is formed.
- **Centrifugation:** 1800-2200g for 8-15 minutes at 18-25° C.



| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|---------|----------------|-----------|----------|-----------|----------|
| GDN050SGC | 5.0 | 13X75 | Clot Activator | Gel | PP | | Push Cap |
| GDN060SGC | 6.0 | 13X100 | Clot Activator | Gel | PP | | Push Cap |
| GDN085SGC | 8.5 | 18X100 | Clot Activator | Gel | PP | | Push Cap |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

Micro Blood Collection Tubes

Used for pediatric & finger-prick blood sample collections. Made from medical grade polypropylene (PP), virtually unbreakable and highly transparent.



Description

- Ultrasonically sprayed additives for different clinical lab tests.
- Precision molded snap cap maintains a tight seal during procedures and allows to open by single hand.
- Caps are available in bright attractive colors for sample identification purposes.
- An integrated scoop on the tube for easier micro blood collection.
- Marked with fill lines, lot numbers and expiration date on labels.
- Holds 250 to 500µL of blood.
- Immediately after draw, gently invert the tube in order to mix with additives.
- To obtain serum with clot activator, tubes should be allowed to stand for 30 minutes or until the clot is formed.

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|----------|-----------|---------|---------------------------------|-----------|----------|-----------|----------|
| GD005A | 0.5 | 10X45 | NONE | / | PP | | Snap Cap |
| GD005CA | 0.5 | 10X45 | Clot Activator | / | PP | | Snap Cap |
| GD005EK2 | 0.5 | 10X45 | EDTA K2 | / | PP | | Snap Cap |
| GD005EK3 | 0.5 | 10X45 | EDTA K3 | / | PP | | Snap Cap |
| GD005LH | 0.5 | 10X45 | Heparin Lithium | / | PP | | Snap Cap |
| GD005SH | 0.5 | 10X45 | Heparin Sodium | / | PP | | Snap Cap |
| GD005SFH | 0.5 | 10X45 | Heparin Sodium+ Sodium Fluoride | / | PP | | Snap Cap |
| GD005SGC | 0.5 | 10X45 | Clot Activator | Gel | PP | | Snap Cap |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.



Micro Blood Collection Tubes



Skirted tube with false bottom to make a full-size blood collection tube, with dimensions as 13 X 75mm; Compatible with automated hematology analyzers without the need for a tube adapter.

Used for pediatric & finger-prick blood sample collections.

Description

- Ultrasonically sprayed additives for different clinical lab tests.
- Precision molded snap cap maintains a tight seal during procedures and allows to open by single hand. They are available in bright attractive colors for sample identification purposes.
- An integrated scoop on the tube for easier micro blood collection.
- Marked with fill lines, lot numbers and expiration date on labels.
- Made from medical grade polypropylene (PP), virtually unbreakable and highly transparent.
- Immediately after draw, gently invert the tube in order to mix with additives.
- To obtain serum with clot activator, tubes should be allowed to stand for 30 minutes or until the clot is formed.

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|------------|-----------|---------|---------------------------------|-----------|----------|-----------|----------|
| GDNM005A | 0.5 | 13X75 | NONE | / | PP | | Snap Cap |
| GDNM005CA | 0.5 | 13X75 | Clot Activator | / | PP | | Snap Cap |
| GDNM005EK2 | 0.5 | 13X75 | EDTA K2 | / | PP | | Snap Cap |
| GDNM005EK3 | 0.5 | 13X75 | EDTA K3 | / | PP | | Snap Cap |
| GDNM005LH | 0.5 | 13X75 | Heparin Lithium | / | PP | | Snap Cap |
| GDNM005SH | 0.5 | 13X75 | Heparin Sodium | / | PP | | Snap Cap |
| GDNM005SFH | 0.5 | 13X75 | Heparin Sodium+ Sodium Fluoride | / | PP | | Snap Cap |
| GDNM005SGC | 0.5 | 13X75 | Clot Activator | Gel | PP | | Snap Cap |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

Capillary Blood Collection Tube

For collecting, transport and processing capillary blood, especially in children, geriatric patients when small amount of blood sample is required.



Description

- Ultrasonically sprayed additives for different clinical lab tests.
- Small volume blood collection use.
- Micro Tube with central crossed rubber cap for inserting capillary tubes, allowing appropriate capillary to be inserted without removing the cap.
- Immediately after draw, gently invert the tube in order to mix with additives.
- To obtain serum with clot activator, tubes should be allowed to stand for 30 minutes or until the clot is formed.

| Cat. No | Volume/ml | Size/mm | Additive | Separator | Material | Cap Color | Closure |
|-----------|-----------|---------|---------------------------------|-----------|----------|-----------|----------|
| GDC005A | 0.5 | 8X42 | NONE | / | PP | | Snap Cap |
| GDC005CA | 0.5 | 8X42 | Clot Activator | / | PP | | Snap Cap |
| GDC005EK2 | 0.5 | 8X42 | EDTA K2 | / | PP | | Snap Cap |
| GDC005EK3 | 0.5 | 8X42 | EDTA K3 | / | PP | | Snap Cap |
| GDC005LH | 0.5 | 8X42 | Heparin Lithium | / | PP | | Snap Cap |
| GDC005SH | 0.5 | 8X42 | Heparin Sodium | / | PP | | Snap Cap |
| GDC005SFH | 0.5 | 8X42 | Heparin Sodium+ Sodium Fluoride | / | PP | | Snap Cap |
| GDC005SGC | 0.5 | 8X42 | Clot Activator | Gel | PP | | Snap Cap |
| GDC005SC | 0.5 | 8X42 | Sodium Citrate | / | PP | | Snap Cap |

- All tubes are supplied in 1000pcs or 1200 pcs/ctn.

VACUUM BLOOD COLLECTION SYSTEM




Butterfly Needles

Safety blood collection needle devices specialized for venipuncture.

Description

- Accessing a superficial vein or artery for intravenous injection or phlebotomy. It consists, from front to rear, of a Hypodermic needle, two bilateral flexible "wings".
- Flexible small-bore transparent tubing (often 20-35cm long).
- Transparent tubes for optimal observation of venous blood return.
- All butterfly needles are supplied in 2000 pcs/ctn.



| Cat. No | Specification | Color |
|----------|---------------|-------|
| GD-BF21G | 21Gx3/4" | Green |
| GD-BF22G | 22Gx3/4" | Black |
| GD-BF23G | 23Gx3/4" | Blue |

Pen Type Blood Collection Needles

Used for blood collection.



Description

- Easily identified by different color .
- Thin-walled design, big inner dia, so the flow rate is high.
- Pyrogen free.
- Non-toxic.

| Cat. No | Size | Color |
|---------|------------|--------|
| GDN18G | 18GX1 1/2" | Pink |
| GDN20G | 20GX1 1/2" | Yellow |
| GDN21G | 21GX1 1/2" | Green |
| GDN21G2 | 21GX1 1/4" | Green |
| GDN22G | 22GX1 1/2" | Black |
| GDN22G | 22GX1 1/4" | Black |
| GDN23G | 23GX1 1/2" | Blue |

- All needles are supplied in 4000 pcs/ctn.

Flashback Blood Collection Needles

Flashback blood collection needles avoid the potential multiple sticks of the same patient with instant confirmation of needle placement, through a unique visual flashback indication.

| Cat. No | Specification | Color |
|----------|--------------------------|-------|
| GDN21GFB | Flash Back Needle 21Gx1" | Green |
| GDN22GFB | Flash Back Needle 22Gx1" | Black |

- All flash back needles are supplied in 5000 pcs/ctn.



Blood Collection Tube Holder

Blood collection tube holder for taking blood samples.

Description

- Compatible with standard, multiple sample blood collection needles.

| Cat. No | Specification | Qty/ctn |
|---------|------------------------------|---------|
| GD-H | Holder, PP | 2000 |
| GD-H1 | Holder with safety cover, PP | 2000 |







Sharps Container




- **Material:** Polypropylene, PVC-free.
- Standard disposable sharps container used to collect sharps waste, including syringes, infusion sets, hypodermic needles, surgical needles, etc.
- **Available in a variety of sizes:**
- 1L, 3L, 5L, 7L, 13L; 1 gallon, 2 gallons, 6 gallons; 1 quart, 2 quarts, 5 quarts.
- The container is made of puncture-resistant medical plastics; the lid is flip-top or rotating for easy collection and transport. FDA 510(k) approval.



E1



X



F

| Item | Usable Capacity | Color | Type | Closure | A1 | B1 | C1 | D1 | E1 Unit(mm) | Qty/ctn |
|-----------|-----------------|--------|---------------|------------|-------|-------|-------|-------|-------------|---------|
| GDC-03-X | 3L | Yellow | Free Standing | Rotary Lid | 184.4 | 183.3 | 141.8 | 154.3 | 184.3 | 48 |
| GDC-05-X | 5L | Yellow | Free Standing | Rotary Lid | 221 | 213.3 | 173.7 | 195.9 | 221 | 32 |
| GDC-07-X | 7L | Yellow | Free Standing | Rotary Lid | 241.2 | 224 | 194 | 220.4 | 241.2 | 24 |
| GDC-013-X | 13L | Yellow | Side Knob | Rotary Lid | 300.2 | 394.3 | 220.7 | 253.6 | 300.2 | 20 |


| Item | Usable Capacity | Color | Type | Closure | A1 | B1 | C1 | D1 | E1 Unit(mm) | Qty/ctn |
|----------|-----------------|--------|---------------|----------|-------|-------|-------|-------|-------------|---------|
| GDC-01-F | 1L | Yellow | Free Standing | Flip Lid | 102.5 | 166.4 | 70.4 | 122 | 102.5 | 60 |
| GDC-03-F | 3L | Yellow | Free Standing | Flip Lid | 158.1 | 190.6 | 112.8 | 188.4 | 58.1 | 48 |
| GDC-05-F | 5L | Yellow | Free Standing | Flip Lid | 199 | 210.1 | 149.1 | 226.2 | 199 | 32 |
| GDC-07-F | 7L | Yellow | Free Standing | Flip Lid | 223.9 | 223.6 | 170.5 | 241.4 | 220.9 | 24 |




Sharps Container




SC-01Q




SC-02Q





SC-05Q



SC-051Q





SC-01G





SCT50



SC-02G




SC-06G



| Item | Usable Capacity | Color | Type | Closure | Qty/ctn |
|---------|-----------------|-------------|-----------------------------------|----------------------|---------|
| SC-01Q | 1 Quart | Red | One-handed | Flid | 60 |
| SC-02Q | 2 Quart | Red | One-handed | Flid Lid | 60 |
| SC-05Q | 5 Quart | Red | Wall Mount | Horizontal Entry Lid | 20 |
| SC-051Q | 5 Quart | Red | Wall Mount | Horizontal Entry Lid | 20 |
| SC-01G | 1 Gallon | Red | Wall Mount | Slide Lid | 24 |
| SC-02G | 2 Gallon | Red | Wall Mount | Slide Lid | 24 |
| SC-06G | 6 Gallon | Red | Side Knob | Slide Lid | 12 |
| SCT50 | 50ml | Transparent | Portable Sharps Shuttle Container | Lock-safe cap | 24 |



SPECIMEN CONTAINERS




Stool specimen collection

Those containers are used for the collection of stool specimens. It provides easy stool collection and closed container transport in a clean, convenient and safe way, effectively preventing the drying, spillage and leakage of stool specimens. It offers easy specimen labeling, bulk transport and sampling during testing.



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--------------------------------------|----------|---------------|-----------------|
| K1011 | 20ml stool container with stick | 1000 | 100pcs*10bags | 48*32*32 |
| K1012 | 25ml stool culture bottle with stick | 650 | 25pcs*26bags | 44*30*35 |
| K1013 | 60ml Stool container with stick | 500 | 1pc*500bags | 49*34*39 |
| K1013-1 | 60ml Stool container without stick | 500 | 1pc*500bags | 49*34*39 |
| H1016-6 | 60ml Stool container with spoon | 500 | 1pc*500bags | 48*32*52 |

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|--------------|-----------------|
| K1012-1 | 25ml stool container | 650 | 25pcs*26bags | 46*30*35 |
| K1014 | 5ml stool collection with spoon | 2000 | 25pcs*80bags | 49*33*33 |
| K1014-1 | 5ml stool collection without spoon | 2000 | 25pcs*80bags | 49*33*33 |
| K1017-2 | 25ml universal container with spoon, container made of PP | 650 | 25pcs*26bags | 56*30*35 |
| K1017-3 | 25ml universal container with spoon, container made of PS | 650 | 25pcs*26bags | 56*30*35 |
| K1017-4 | 25ml universal container without spoon, container made of PP | 650 | 25pcs*26bags | 56*30*35 |
| K1018-2 | 20ml stool container with spoon | 1500 | in bulk | 42.5*38*42 |
| K1018-3 | 20ml stool container with spoon | 650 | in bulk | 37.5*29.5*33 |



Collection of sputum cup specimens


This container is used for the collection of sputum specimens. It provides easy sputum collection and closed container transport in a clean, convenient and safe way.



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|---------------|-----------------|
| K1015 | 30ml sputum container | 500 | 1pc*500bags | 48*32*34 |
| K1016 | 40ml sputum container | 500 | 5pcs*100bags | 39*39*39 |
| K1017 | 25ml universal container | 650 | 25pcs*26bags | 46*30*32 |
| K1017-1 | 25ml sputum container (skirted bottom) | 650 | 25pcs*26bags | 50*30*35 |
| K1018 | 30ml sputum container | 1000 | 10pcs*100bags | 56.5*31*57 |




| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|-----------------------|----------|----------------|-----------------|
| K1018-1 | 20ml sputum container | 1500 | in bulk | 42.5*38*42.5 |
| K1019 | Sputum cassette | 1000 | 100pcs*10bags | 48*32*31 |
| K1023 | Mucus trap, 40CC | 50 | 1pc/bag*50bags | 40*27*20 |
| K1024 | Mucus trap, 80CC | 50 | 1pc/bag*50bags | 40*27*20 |




Stool specimen collection


This container is used for the collection of urine specimens. It is made of medical grade plastics with a user-friendly design. It has an evenly thick wall that enables easy urine collection and closed transport and testing. It thus facilitates the specimen transport for testing and its use with testing instruments, and prevents contamination. (The product use is clean, convenient and safe, which provides a guarantee for the normalization and standardization of urine analysis from the source.)



K1025




K1010




K1009




H1015



H1016



H1016-1



H1016-2



H1016-5



H1021




H1025



K1026




K1027



GD-H3



GD-H4



GD080UAC GD0100UA GD0100UB GD080UBC GD0100C&S



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|---------------|-----------------|
| K1009 | Soft urine cup white | 10000 | 500pcs*20bags | 47*27*33 |
| K1010 | Soft urine cup transparent | 10000 | 500pcs*20bags | 47*27*33 |
| K1025 | 30ml urine container with mouth on cap | 500 | 1pc*500bags | 48*32*34 |
| H1015 | 100ml-120ml specimen container Individual packing | 500 | 1pc*500bags | 54*48.5*61 |
| H1016 | 60ml specimen container | 500 | 1pc*500bags | 48*32*52 |
| H1016-1 | 60ml specimen container, container made of PS | 500 | 1pc*500bags | 43*43*33 |
| H1016-2 | 90ml specimen container | 500 | 1pc*500bags | 53*46*44 |
| H1016-5 | 40ml renovate urine container | 500 | in bulk | 48*31.5*52 |

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|-----------|---|----------|-------------------|-----------------|
| H1021 | 140ml specimen container | 500 | 1pc*500bags | 58*58*57 |
| H1025 | Specimen Container,80CC | 500 | 25pcs*20bags | 69*26*33 |
| K1026 | Urine transfer device for urine tube, 100ml-120ml | 250 | 1pc*250bags | 60*31*44 |
| K1027 | Urine transfer device for urine tube, 100ml-120ml | 250 | 1pc*250bags | 60*31*44 |
| GD080UAC | Vacuum Urine Tube, Conical, Rubber stopper, no Additive, 16*100mm, 8ml | 1200 | 100pcs*12racks | 51*40*24 |
| GD0100UA | Vacuum Urine Tube, Yellow cap, no Additive, 16*100mm, 10ml | 1200 | 100pcs*12racks | 51*40*24 |
| GD0100UB | Vacuum Urine Tube, Yellow cap, Boric Acid, 16*100mm, 10ml | 1200 | 100pcs*12racks | 51*40*24 |
| GD080UBC | Vacuum Urine Tube, Conical, Yellow cap, Boric Acid, 16*100mm, 8ml | 1200 | 100pcs*12racks | 51*40*24 |
| GD0100C&S | Vacuum Urine Tube, Yellow cap, Boric Acid Sodium Borate/Formate, 16*100mm, 10ml | 1200 | 100pcs*12racks | 51*40*24 |
| GD-H3 | Urine Transfer Straw, Yellow, PP | 5000 | 50pcs/bag*100bags | 38*35*17 |
| GD-H4 | Urine Transfer Straw, White, PE | 5000 | 50pcs/bag*100bags | 38*35*17 |



95KPA Compliant Specimen Container

Product is specifically designed for the safe, secure collection and transportation of biological specimens for laboratory evaluation



GD5343 GD4834 GD4833 GD4825 GD4815 GD3510

Urine Storage Container for 24 hours

Using the latest nano antibacterial release technology, made of new antibacterial inert plastic, can significantly inhibit the growth of bacteria within 48 hours. It is especially suitable for containers that need to collect urine continuously for 24 hours for analysis of various proteins, hormones, amino acids and other items.




H1026 H1027 H2001 H2002 H1016-3

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|---------------|-----------------|
| GD3510 | 20ml Container with 35mm Cap | 1000 | EA/Case | 40*40*34 |
| GD4815 | 40ml Container with 48mm Cap | 500 | EA/Case | 49*49*39 |
| GD4825 | 60ml Container with 48mm Cap | 500 | EA/Case | 49*49*39.5 |
| GD4833 | 90ml Container with 53mm Cap | 300 | EA/Case | 56*56*29 |
| GD4834 | 90ml Container with 48mm Cap | 400 | EA/Case | 49*49*37.5 |
| GD5343 | 120ml Container with 53mm Cap | 300 | EA/Case | 56*56*32 |
| H1016-3 | 40ml specimen container | 500 | 1pc*500bags | 42*33*42.5 |
| H1026 | Urine storage container for 24 hours, 3000ml | 50 | in bulk | 77*31*51 |
| H1027 | Urine tube, 16*85mm, 10ml | 2000 | 200pcs*10bags | 48*32*39 |
| H2001 | 53mm Metal Lid Specimen Container, Sterile, 100ml | 400 | in bulk | 52x42x48 |
| H2002 | 70mm Metal Lid Specimen Container, Sterile, 100ml | 400 | in bulk | 70x35x50 |



...raw material, its smoothness and transparent finish fully achieve the... standard of similar product. It is produced Class 100,000 cleanroom... no heat source and sterilized by ethylene oxide. Adopt the... easy tear package, which is simple and convenient and... opportunities in the use of products. Each batch of products... observation test of cell growth rate, test of cell... and sterility test to ensure product quality

PETRI DISHES



K1001



K1002



K1003



K1009



K1004-J-3 Vents

K1004-M

K1004-T



K1005

K1006

K1007

K1008



K1004-J



K1020

K1021

K1022



K1032-L

K1032-S



K1033

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|----------------|------------------------------------|----------|-------------------|-----------------|
| K1001 | 150*15mm Petri Dish | 180 | 10pcs/bag*18bags | 45*35*45 |
| K1002 | 90*20mm Petri Dish | 500 | 10pcs/bag*50bags | 47*47*44 |
| K1004-J | 90*15mm no vents Petri Dish | 500 | 10pcs/bag*50bags | 47*47*34 |
| K1004-J-3Vents | 90*15mm 3-vents Petri Dish | 500 | 10pcs/bag*50bags | 47*47*34 |
| K1004-M | 90*15mm for machine use Petri Dish | 500 | 10pcs/bag*50bags | 47*47*34 |
| K1004-T | 90*15mm 3-vents Petri Dish | 500 | 10pcs/bag*50bags | 47*47*34 |
| K1020 | Inoculation Loops, 10ul, Blue | 3000 | 10pcs/bag*300bags | 38*28*20 |
| K1021 | Inoculation Loops, 1ul, White | 3000 | 10pcs/bag*300bags | 38*28*20 |
| K1022 | Inoculation Loops, Needles, Red | 3000 | 10pcs/bag*300bags | 38*28*20 |




| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|-------------------------------------|----------|--------------------|-----------------|
| K1003 | 90*15mm Petri Dish Two-Compartments | 500 | 10pcs/bag*50bags | 47*34*47 |
| K1005 | 70*15mm Petri Dish | 600 | 10pcs/bag*60bags | 44*27*35 |
| K1006 | 65*15mm Petri Dish | 650 | 10pcs/bag*65bags | 49*34*28 |
| K1007 | 55*15mm Petri Dish | 1000 | 10pcs/bag*100bags | 46*30*35 |
| K1008 | 30*15mm Petri Dish | 2000 | 10pcs/bag*200bags | 49*33*35 |
| K1009 | 14cm Petri Dish | 180 | 10pcs/bag*18bags | 45*39*45 |
| K1032-S | 25cm Cell Scaper | 100 | Individual Packing | 33*27*11 |
| K1032-L | 39cm Cell Scaper | 100 | Individual Packing | 53*27*11 |
| K1033 | 18cm Cell Lifter | 100 | Individual Packing | 31*20*11 |



NURSING PLASTICWARE SERIES

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|-------------------------------|-------------------|
| P1001-L | Vaginal Speculum, Large | 100 | 10pcs/bag*10bags | 49*24*38 |
| P1001-M | Vaginal Speculum Middle | 100 | 10pcs/bag*10bags | 38*23*39 |
| P1001-S | Vaginal Speculum, Small | 100 | 10pcs/bag*10bags | 38*23*39 |
| P1002-L | Vaginal Speculum with light source , Large | 80 | 8pcs/bag*10bags | 56*40*30 |
| P1002-M | Vaginal Speculum with light source, Middle | 100 | 10pcs/bag*10bags | 56*40*30 |
| P1002-S | Vaginal Speculum with light source, Small | 100 | 10pcs/bag*10bags | 56*40*30 |
| P1003 | Transparent dressing bowl(Sterile) | 1500 | 2pcs/packet*50pcs/bags*15bags | 48*32*32 |
| P1004-1 | Dressing basin(soft) | 1000 | 2pcs/packet*50pcs/bags*10bags | 49*33*33 |
| P1004-2 | Dressing basin(hard) | 80 | 2pcs/bag*40bags | 40*30*25 |
| P1008 | Tongue Depressor | 2800 | 200pcs/bag*14bags | 50*35*37 |
| P1010 | Mouth piece for Endoscope(Large/Medium) | 1000 | 100pcs/bag*10bags | 55*31*43/55*31*38 |


NURSING PLASTICWARE SERIES



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------------|----------------------------------|----------|-------------------|-----------------|
| P1012 a/b/c/e | Tweezers 12.5cm Long | 2500 | 250pcs/bag*10bags | 48*32*32 |
| P1012-1 | Fine-tipped Tweezers 10cm Long | 2500 | 250pcs/bag*10bags | 48*32*44 |
| P1013 | Tweezers with tooth 12cm Long | 2500 | 250pcs/bag*10bags | 48*32*32 |
| P1014 | Long Tweezers 13cm Long | 2000 | 200pcs/bag*10bags | 48*32*32 |
| P1015 | Tweezers without tooth 12cm Long | 2500 | 250pcs/bag*10bags | 48*32*32 |
| P1016 | Blue Tweezers 13cm Long | 2000 | 200pcs/bag*10bags | 48*32*44 |
| P1020 | Sponge Holder 23cm Long | 500 | 50pcs/bag*10bags | 52*38*32 |
| P1021 | Sponge Holder 18cm Long | 500 | 50pcs/bag*10bags | 42*38*32 |
| P1038 | Blue Tweezers 13cm Long, ABS | 2000 | 200pcs/bag*10bags | 48*32*44 |



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|------------------------------------|----------|-------------------|-----------------|
| P1022 | Sponge Holder 11cm Long | 1200 | 100pcs/bag*12bags | 49*33*33 |
| P1023 | Sponge Holder 10cm Long | 1200 | 100pcs/bag*12bags | 49*33*33 |
| P1024 | Sponge Holder 14.3cm Long | 1200 | 100pcs/bag*12bags | 48*32*42 |
| P1025 | Sponge Holder 13.9cm Long | 1200 | 100pcs/bag*12bags | 48*30*39 |
| P1033 | Scissors, Blue, ABS | 1000 | 100pcs/bag*10bags | 49*33*33 |
| P1004 | Scissors, Blue | 1200 | 100pcs/bag*12bags | 49*33*33 |
| P1026 | Disposable umbilical cord scissors | | | |

NURSING PLASTICWARE SERIES



CENTRIFUGE TUBE CRYO VIALS SERIES

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|-------------------------------------|----------|-------------------------|-----------------|
| K1034 | 2.5mm Ear Speculum | 8500 | 34pcs/bag*25bags*10bags | 40*26*27 |
| K1035 | 4.25mm Ear Speculum | 8500 | 34pcs/bag*25bags*10bags | 40*26*27 |
| P1030 | Anus speculum, Large, Middle, Small | 100 | 25pcs/bag*4bags | 48*30*32 |
| P1039 | Anus speculum | 100 | 25pcs/bag*4bags | 30*22*27 |
| K1028 | Ear curettes,White | 500 | 1pc*50pcs/box*10boxes | 31*18.5*13.5 |
| K1029 | Ear curettes,Green | 500 | 1pc*50pcs/box*10boxes | 31*18.5*13.5 |
| K1030 | Ear curettes,Pink | 500 | 1pc*50pcs/box*10boxes | 31*18.5*13.5 |
| K1031 | Ear curettes,Grey | 500 | 1pc*50pcs/box*10boxes | 31*18.5*13.5 |
| K1036 | Ear curettes,Blue | 500 | 1pc*50pcs/box*10boxes | 31*18.5*13.5 |
| P1045 | Light handle cover ,green | 1000 | in bulk | 58*40*45 |
| P1046 | Light handle cover ,Purple | 250 | in bulk | 56*36*47 |



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|----------|--|----------|---------------|-----------------|
| D1001 | 50ml conical with graduation | 500 | 50pcs*10bags | 48*32*43 |
| D1002 | 50ml self-standing | 500 | 50pcs*10bags | 70*40*25 |
| D1003 | 15ml conical with graduation | 1000 | 50pcs*20bags | 46*29*36 |
| D1003-PS | 15ml conical with graduation, container made of PS | 1000 | 50pcs*20bags | 57*27*27 |
| D2003 | 15ml conical with graduation | 1000 | 50pcs*20bags | 46*29*36 |
| D1004 | 50ml | 500 | 50pcs*10bags | 49*33*33 |
| D1005 | 15ml | 1500 | 150pcs*10bags | 49*33*33 |
| D1006 | 10ml | 2000 | 200pcs*10bags | 49*33*33 |
| D1008 | 5ml with cap | 3000 | 300pcs*10bags | 48*32*42 |

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|----------------------------|----------|----------------|-----------------|
| D1009 | 2ml | 5000 | 500pcs*10bags | 48*32*32 |
| D1010 | 1.5ml with transparent cap | 6000 | 500pcs*12bags | 49*33*33 |
| D1011 | 0.5ml | 12000 | 1000pcs*12bags | 49*33*33 |
| D1012 | 0.2ml | 30000 | 1000pcs*30bags | 48*32*32 |
| D1017 | 1.5 ml linked green cap | 6000 | 500pcs*12bags | 49*33*33 |
| D1018 | 1.5ml linked red cap | 6000 | 500pcs*12bags | 49*33*33 |
| D1019 | 1.5ml blue cap | 6000 | 1000pcs*6bags | 49*33*33 |
| D1020 | 1.5 ml red cap | 6000 | 1000pcs*6bags | 49*33*33 |
| D1039 | 2ml conical | 5000 | 500pcs*10bags | 48*32*32 |



D1013

D1015

D1013/D1015

- Material: Medical grade polypropylene (PP).
- Withstands temperatures from -20°C to 100°C; can be used for liquid nitrogen (-196°C) freezing in static conditions.
- Sterilization: E-Beam sterilization or EO sterilization or non-sterile.



D1015-20

D1015-15

D1015-05

D1015-10

D3027

D3028

D3029

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|----------|----------------------|----------|---------------|-----------------|
| D1013 | 5ml+cap | 2000 | 200pcs*10bags | 49*33*33 |
| D1015 | 1.8ml | 5000 | 500pcs*10bags | 48*32*32 |
| D1015-20 | Cryo vial 2ml | 5000 | 500pcs*10bags | 48*32*32 |
| D1015-15 | Cryo vial 1.5ml | 5000 | 500pcs*10bags | 48*32*32 |
| D1015-05 | Cryo vial 0.5ml | 5000 | 500pcs*10bags | 48*32*32 |
| D1015-10 | Cryo vial 1ml | 5000 | 500pcs*10bags | 48*32*32 |
| D3027 | Cryo vial 0.5ml | 5000 | 500pcs*10bags | 48*32*32 |
| D3028 | Cryo vial 1.5ml | 5000 | 500pcs*10bags | 48*32*32 |
| D3029 | Cryo vial 2.0ml | 5000 | 500pcs*10bags | 48*32*32 |



D1038

D1037

D3056

D3033

D3034

D3035

D3036

D3037

D3038

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|--|-----------------|
| D3033 | Cryo vial 0.5ml | 5000 | 500pcs*10bags | |
| D3034 | Cryo vial 1.0ml | 5000 | 500pcs*10bags | |
| D3035 | Cryo vial 1.8ml | 5000 | 500pcs*10bags | |
| D3036 | Cryo vial 1.5ml | 5000 | 500pcs*10bags | |
| D3037 | Cryo vial 1.8ml | 5000 | 500pcs*10bags | |
| D3038 | Cryo vial 5.0ml | 5000 | 500pcs*10bags | |
| D1037 | 5ml transport vials | 2000 | Caps:1000pcs/bag*2bags Tubes:200pcs/bag*10bags | 48*48*42 |
| D1038 | 10ml transport vials | 1600 | Caps:800pcs/bag*2bags Tubes:200pcs/bag*8bags | 47*42*34 |
| D3056 | 0.5ml cryovial,3codes (Barcode,digital code,and 2D) | 1920 | 96pcs/box*20boxes | |
| D3057 | 0.5ml cryovial,3codes (Barcode,digital code,and 2D) | 2000 | 100pcs/bag*20bags | |




PCR Eight connected pipe

- Material: Medical grade polypropylene.
- The tube is transparent for easy sample observation.
- Dished cap: ensures a tight fit to effectively prevent the evaporation and loss of samples;
- Flat cap: ideal for quantitative real-time PCR experiments. The tube cap is specially lengthened for inner seal and the cap body is sealed, effectively preventing accidental opening.
- Customized support: Can be DNA/RNA-free, pyrogen-free and endotoxin-free.
- Operating temperature: -20°C to 100°C.
- Fit for common brand PCR instruments; withstand centrifugal force of 6,000G in a fully supported rotor.



D1028



D1029





| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|------------------------------|-----------------|
| D1028 | 8-Strip tubes for PCR 0.1ml milk white with cap | 1250 | 125pcs/bag/small box*10boxes | 47*30*15 |
| D1029 | 8-Strip tubes for PCR 0.2ml transparent with cap | 1250 | 125pcs/bag/small box*10boxes | 47*30*15 |

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|------------|-----------------|
| D1021 | 98 wells centrifugation tube rack for 0.5ml | 40 | 1pc*40bags | 48*33*32 |
| D1022 | 72 wells centrifugation tube rack for 1.5ml | 40 | 1pc*40bags | 48*33*32 |
| D1023 | 81 wells centrifugation tube rack for 1.5ml | 40 | 1pc*40bags | 48*33*32 |
| D1024 | Two purpose tube rack | 40 | 1pc*40bags | 48*33*32 |
| D1027 | Two-purpose centrifuge tube rack(front 1.5ml, reverse side 0.5ml) | 50 | 1pc*50bags | 48*33*32 |
| H1013 | Centrifugation rack for 50ml, 15ml | 50 | 1pc*50bags | 44*28*37 |







PIPETTE TIPS

TO ENSURE the accuracy and precision of the experiment. The precision pipette tip provides optimum reproducibility and accuracy. Standard tip and tip with filter can be disinfected and sterilized in high pressure. The unique filter element silicified tip does not hang the residual liquid, which guarantees the accuracy of liquid transfer to the utmost. It almost meets and fits the requirements of most commercially available pipettes. It is also equipped with a variety of tip boxes, in high temperature and pressure. Match with a variety of centrifuge tube holders.



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|--------------------|-----------------|
| E1010 | 10µl φ5*32mm Fit For Gilson | 60000 | 1000pcs/bag*60bags | 49*33*33 |
| E1011 | 10µl-200µl φ6*50mm (Yellow) Fit For Gilson | 30000 | 1000pcs/bag*30bags | 48*32*32 |
| E1012 | 10µl-200µl φ6*50mm (Graduation) Fit For Gilson | 30000 | 1000pcs/bag*30bags | 48*32*32 |
| E1013 | 200µl-1000µl φ8*71mm Fit For Gilson | 10000 | 500pcs/bag*20bags | 48*32*40 |
| E1014 | 10µl φ5*32mm Fit For Finland | 60000 | 1000pcs/bag*60bags | 49*33*33 |
| E1015 | 10µl-200µl φ6*50mm Fit For Finland | 30000 | 1000pcs/bag*30bags | 48*32*32 |
| E1016 | 10µl-250µl φ6*52mm Fit For Finland | 25000 | 1000pcs/bag*25bags | 48*32*32 |
| E1017 | 200µl-1000µl φ8*71mm Fit For Finland | 10000 | 500pcs/bag*20bags | 48*32*40 |

PIPETTE TIPS



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|--------------------|-----------------|
| E1007 | 200µl-1000µl Fit For Eppendorf | 10000 | 500pcs/bag*20bags | 48*32*40 |
| E1008 | 10µl-300µl Fit For Eppendorf | 30000 | 1000pcs/bag*30bags | 48*32*44 |
| E1009 | 10µl φ4*49mm Fit For Eppendorf | 60000 | 1000pcs/bag*60bags | 48*32*32 |
| E1018 | 200µl-1000µl φ8*71mm Fit For Switzerland | 10000 | 500pcs/bag*20bags | 48*32*40 |
| E1019 | 10µl-200µl φ6*50mm Fit For Switzerland | 30000 | 1000pcs/bag*30bags | 48*32*32 |
| E1010-1 | 10µl φ5*32mm With Filter Fit For Gilson | 60000 | 1000pcs/bag*60bags | 49*33*33 |
| E1020-1 | 10µl-1000µl With Filter Fit For Gilson | 10000 | 500pcs/bag*20bags | 49*33*40 |
| E1021-1 | 10µl-200µl φ6*50mm With Filter Fit For Gilson | 30000 | 1000pcs/bag*30bags | 48*32*40 |
| E1024 | Tips Box 96wells For 200µl | 68 | in bulk | 48*32*32 |
| E1025 | Tips Box 96wells For 200µl&10µl | 68 | in bulk | 48*32*32 |
| E1027 | Fit for Uranus Ae120 | 9600 | 192pcs/box*50boxes | 66.5*46.5*28 |

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|------------------------------|-----------------|
| E1032 | Fit For ADDCON | 10000 | 100pcs/bag*100bags | 57*35*22 |
| E1033 | Fit For Siemens | 8640 | 24box/carton | 60*35*38 |
| E1034 | Fit For TECAN Freedom EVO 200 MAC 96 Workstation | 24000 | (1000pcs*2bags) /box*12boxes | 56*52*35 |
| E1035 | TECAN Freedom EVO MCA Workstaion | 24000 | (1000pcs*2bags) /box*12boxes | 56*52*35 |
| E4503 | 50-300µl Fit For Finland | 11520 | 384pcs/box*30boxs | 73*36*38 |
| E3100 | 1000µl Fit For TECAN | 4608 | 192pcs/box*24boxs | 57*46*45.5 |
| E3101 | 200µl Fit For TECAN | 4608 | 192pcs/box*24boxs | 57*46*32 |

Pyrogen-free Tips

- Material: Polypropylene (PP).
- These tips have a hydrophobic filter to prevent cross-contamination of aerosol and liquid. The low-adsorption tips are suitable for handling viscous liquids and liquids with low surface tension, greatly reducing the liquid retention on the wall to ensure the accuracy of experimental results. They can be directly used in genetic testing.



E1010-1D



E1021-1D



E1020-1D



E1020-1D1250




E1022

E1023




CULTURE PLATE & PIPETTE SERIES

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|--------------|---|----------|---|-----------------|
| E1010-1D | Universal tips with filter for dispensers 10 µl | 4800 | 96pcs/plastic box*10boxes/carboard box*5carboard box/Case | 48*27*32 |
| E1020-1D | Universal tips with filter for dispensers 1000 µl | 4800 | 96pcs/plastic box*10boxes/carboard box*5carboard box/Case | 48.5*27*61.5 |
| E1020-1D1250 | Universal tips with filter for dispensers 1250 µl | 4800 | 96pcs/plastic box*10boxes/carboard box*5carboard box/Case | 48.5*27*61.5 |
| E1021-1D100 | Universal tips with filter for dispensers 100 µl | 4800 | 96pcs/plastic box*10boxes/carboard box*5carboard box/Case | 48*27*36 |
| E1021-1D | Universal tips with filter for dispensers 200 µl | 4800 | 96pcs/plastic box*10boxes/carboard box*5carboard box/Case | 48*27*36 |
| E1022 | Tips Box 100wells For 1000µl | 30 | in bulk | 48*32*32 |
| E1023 | Tips Box 200wells For 200µl | 40 | in bulk | 48*32*32 |



96 Wells Plate

- Specially designed with 96 wells; each well has a capacity of 0.2ml.
- Precisely molded for uniformly thin walls that allows well-to-well consistency and fast heat transfer.
- This plate has a U-shaped bottom or a V-shaped bottom. It features a flexible design with removable boards that allow a combination of 12/24/48/72/96 wells as needed.
- Surface marked with letters and numbers for easy positioning and identification.



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|-------------|-----------------|
| F1003-U | 96 wells Plate (U bottom) | 270 | 1pc*270bags | 48*32*37 |
| F1003-F | 96 wells Plate (Flat bottom) | 270 | 1pc*270bags | 48*32*37 |
| F1003-V | 96 wells Plate (V bottom) | 270 | 1pc*270bags | 48*32*37 |
| F1004 | 24 wells Plate | 170 | 1pc*170bags | 43*27*21 |
| F1008 | 20 Wells plate | 180 | 1pc*180bags | 49*33*33 |
| F1021 | 96 wells detachable plate, suitable for 8 or 12 wells strip | 270 | 1pc*270bags | 43*27*21 |

Square Media Bottles

- Flask material: PC or PET; cap material: polyethylene
- Clear graduation marks with an error of ± 5
- Features an anti-reverse flow mouth; available in five sizes: 125ml, 150ml, 250ml, 500ml, 1000ml
- Sterilization: E-Beam sterilization or EO sterilization; D-Nase free, R-Nase free, pyrogen free



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|----------|--|----------|--------------------|-----------------|
| PET250W | PET Square Media Bottles 250ml, White Cap | 100 | Individual Package | 68*34*32 |
| PET500W | PET Square Media Bottles 500ml, White Cap | 100 | Individual Package | 79*40*38 |
| PET1000W | PET Square Media Bottles 1000ml, White Cap | 50 | Individual Package | 47*47*46 |



Embedding Cassette

Embedding cassette is made of durable and porous molded plastic, which can secure the samples stored in liquid and completely avoid the chemical reaction of histological solvent and provide a variety of different colors for easily identifying different tissue samples.




- Material: Engineering plastic POM with good chemical resistance to all histological solvents.
- Writing area is at 45° available for inkjet printing and laser printing.
- Available in white, pink, green, yellow, orange, and blue.

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|----------------------------------|----------|------------------|-----------------|
| F1009 | Embedding Cassette | 5000 | 500pcs*10bags | 48*32*48 |
| F1009H | Embedding Cassette in Tube Racks | 3000 | 75ea/tube*40tube | |
| F1010 | Embedding Cassette White | 4000 | 500pcs*8bags | 48*32*32 |
| F1011 | Embedding Cassette Green | 2700 | 300pcs*9bags | 48*32*32 |
| F1012 | Embedding Cassette Hinged Lid | 5000 | 500pcs*10bags | |


Transfer Pipettes

- Plastic transfer pipettes make up for the fragility of glass pipettes, ending the history of transferring high-risk contagious liquids by glass pipettes and rubber plugs.
- Smooth inner wall reduces the damage of tissue structure and cells in the liquid, and allows quick, safe, and quantitative transfer of the liquid specimen.
- Cold-resistant; can be sterilized at high temperatures.



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|-------------------------------|----------|------------------------|-----------------|
| F1014 | 25 µl & 35µl Transfer pipette | 16000 | 2000pcs*8bags | 48*32*32 |
| F1015 | 1ml Transfer pipette | 5000 | 250pcs*20bags | 54*34*40 |
| F1016 | 3ml Transfer pipette | 5000 | 250pcs*20bags | 54*34*42 |
| F1016-l | 3ml Transfer pipette, sterile | 2500 | 1pc*100pcs/bags*25bags | 61*41*41 |
| F1016-1 | 1.5ml Transfer pipette | 5000 | 250pcs*20bags | 54*34*38 |
| F1016-2 | 1ml Transfer pipette | 5000 | 250pcs*20bags | 54*34*38 |
| F1016-3 | 3.5ml Transfer pipette | 4000 | 250pcs*16bags | 54*34*30 |

**CULTURE PLATE
PIPETTE SERIES**








| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|----------|------------------------|----------|---------------|-----------------|
| F1016-4 | 1ml Transfer pipette | 5000 | 250pcs*20bags | 54*34*38 |
| F1016-5 | 3.5ml Transfer pipette | 4000 | 250pcs*16bags | 54*34*40 |
| F1016-6 | 2ml Transfer pipette | 5000 | 250pcs*20bags | 54*34*42 |
| F1016-7 | 5ml Transfer pipette | 4000 | 250pcs*16bags | 54*34*40 |
| F1016-8 | 1ml Transfer pipette | 5000 | 250pcs*20bags | 54*34*38 |
| F1016-9 | 2ml Transfer pipette | 5000 | 250pcs*20bags | 54*34*42 |
| F1016-11 | Spherical pipette | 4000 | 250pcs*16bags | 49*33*43 |



TEST TUBES & RACKS SERIES



**TEST TUBES
RACKS SERIES**



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|---------------|-----------------|
| G1030 | Test Tube 16*100mm | 1600 | 200pcs*8bags | 42*33*35.5 |
| G1001 | Test Tube 15*100mm, 10ml | 2000 | 200pcs*10bags | 42*33*35.5 |
| G1002 | Test Tube 13*100mm, 7ml | 3000 | 300pcs*10bags | 43*35*38 |
| G1003 | Test Tube 13*78mm, 5.5ml | 4008 | 334pcs*12bags | 48.5*32.5*31.5 |
| G1004 | Test Tube 12*60mm, 3.5ml | 6000 | 400pcs*15bags | 49.5*33*33 |
| G1005 | Test Tube 12*75mm, 5ml | 4800 | 400pcs*12bags | 49.5*33*33 |
| G1006 | Test Tube 12*78mm, 5.5ml | 4800 | 400pcs*12bags | 49.5*33*33 |
| G1009 | Test Tube 13*100mm, 7ml with cap, with 5ml graduation | 3000 | 300pcs*10bags | 43*33*42 |
| G1010 | Test Tube 12*75mm, 5ml with cap, with 2ml graduation | 4800 | 400pcs*12bags | 49*33*38 |
| G1011 | Test Tube 12*60mm, 3.5ml, with 2ml graduation | 6000 | 400pcs*15bags | 49*33*33 |

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|---------------|-----------------|
| G1012 | Test Tube 12*75mm, 5ml, used for BD flow cytometry | 4008 | 334pcs*12bags | 48*32*32.5 |
| G1013 | Test Tube 12*75mm, 5ml, hexagonal bottom | 4800 | 400pcs*12bags | 48*32*32.5 |
| G1019 | Test Tube 15*50mm, 4.5ml | 3600 | 400pcs*9bags | 48*33*33 |
| G1020 | Test Tube 13*78mm, 5.5ml | 4800 | 400pcs*12bags | 49*32*39 |
| G1021 | Test Tube 12*75mm, 5ml, container made of PS/PP | 4800 | 400pcs*12bags | 49*32*39 |
| G11021 | Test Tube 16*75mm (round bottom) | 2500 | 25pcs*100bags | 49*33*33 |
| G1025 | Test Tube 15*60mm (flat bottom) | 2800 | 5pcs*560bags | 49*33*33 |
| G1026 | Test Tube 16*60mm (round bottom) | 2500 | 25pcs*100bags | 49*33*33 |
| G1027 | Test Tube 15*150mm, 15ml | 1500 | 100pcs*15bags | 49*33*33 |
| G1007 | Test Tube 16*53mm, 6ml with cap | 3000 | 300pcs*10bags | 49.5*33*37 |

**TEST TUBES
RACKS SERIES**




Product Description

- Tube body material: high transparent PS, Cap: PP + silicone film.
- Maximum capacity: 11ml, can withstand 3400 rpm centrifugation.

Use Method

- After collecting urine with a disposable urine cup, the tested person pours the urine into the urine sediment test tube, puts the cap on and sends it to the laboratory, and puts it into the special test tube rack of the urine sediment analyzer. The cap should be unscrewed before starting the test.



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|----------------|-----------------|
| G1016 | Test Tube 12ml with cap | 2800 | 200pcs*14bags | 48*32*52 |
| G1017 | Test Tube 12ml urine cup | 1500 | 150pcs*10bags | 48*33*33 |
| G1017-1 | 12ml urine cup with cap | 1500 | 150pcs*10bags | 48*33*42 |
| G1017-2 | US urine test tube | 1200 | 200pcs*6bags | 46*32*29 |
| G1017-3 | Urine test tube(new type) with cap | 2000 | 100pcs*20bags | 43*35*33 |
| G1017-4 | 12ml urine cup with cap | 1500 | 150pcs*10bags | |
| H1010 | Tube stopper For Φ13mm | 20000 | 2000pcs*10bags | 49*32*42 |
| H1011-1 | Tube stopper For Φ13mm | 18000 | 2000pcs*9bags | 49*32*42 |
| H1012 | Tube stopper For Φ15mm | 15000 | 1000pcs*15bags | 49*32*42 |
| H1014-1 | Tube stopper(universal type),For Φ12mm | 22000 | 2000pcs*11bags | 49*33*33 |
| H1014-5 | Tube stopper(universal type),For Φ12mm—Φ16mm | 15000 | 1000pcs*15bags | 49*32*42 |

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|-----------|---|----------|--------------------------|-----------------|
| G1017-14 | 6ml urine test tube | 2000 | 200pcs/bag*10bags | |
| G1022 | 13*150mm, tube with Swab | 1800 | 150pcs*12bags | 57*35*27 |
| G1023 | 12*75mm, tube with Swab | 2200 | 2pcs/bags*100bags*11bags | 48*32*32 |
| G1026-6cc | Test tube with graduation with cap | 500 | in bulk | 35*29*27 |
| G1026-8cc | Test tube with graduation,14.7X101mm with cap | 500 | in bulk | 31*31*35 |
| G1031 | 16*95mm test tube with pipette | 1500 | 150pcs*10bags | 49*33*33 |
| G1035 | Test Tube 12*75mm,container made of PP | 5000 | 500pcs/bag*10bags | 48*31*38 |

**TEST TUBES
RACKS SERIES**



G10261

G10262

G10263

G10264

G10265

G1029FT

G1029FN

Cary-Blair Transport medium

- It is used to collect and transport specimen of Escherichia coli, Klebsiella, Salmonella and Shigella.
- Specimens can be stocked at 4°C for 48 hours.



Amies and Amies with Charcoal Transport medium

- It is used to collect and transport specimen of Escherichia coli.
- Specimens can be stocked at 4°C for 48 hours.

Stuart and Stuart with Charcoal Transport medium




- It is used to collect and transport specimen of gonococcus.
- Specimens can be stocked at 4°C for 48 hours.

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|--|-----------------|
| G10261 | Transport swab with medium (Amies) | 600 | 50pcs*12bags | 60*40*24 |
| G10262 | Transport swab with medium (Amies with charcoal) | 600 | 50pcs*12bags | 60*40*24 |
| G10263 | Transport swab with medium (Cary Blair) | 600 | 50pcs*12bags | 60*40*24 |
| G10264 | Transport swab with medium (Stuart) | 600 | 50pcs*12bags | 60*40*24 |
| G10265 | Transport swab with medium (Stuart with charcoal) | 600 | 50pcs*12bags | 60*40*24 |
| G1029FT | Oropharyngeal swab with flocking tip for throat sampling,EO sterile | 4000 | 1pc/peel pouch*20pcs/small box*20boxes | 44x36x43 |
| G1029FN | Nasopharyngeal swab with flocking tip for nasal sampling,EO sterile | 4000 | 1pc/peel pouch*20pcs/small box*20boxes | 44x36x43 |



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---------------------------------|----------|-----------------|-----------------|
| H1001 | 40 Wells φ18mm Tube Rack | 40 | 2pcs/bag*20bags | 38*25*25 |
| H1002 | 50Wells φ15mm Tube Rack | 40 | 2pcs/bag*20bags | 38*25*25 |
| H1003 | 50Wells φ13mm Tube Rack | 40 | 2pcs/bag*20bags | 38*25*25 |
| H1004 | Pillar-like φ5*16 Tube Rack | 54 | 1pc*54bags | 48*32*32 |
| H1005 | Pillar-like φ5*10 Tube Rack | 51 | 1pc*51bags | 48*32*32 |
| H1018 | Multi-Function Rack | 50 | 1pc*50bags | 61*41*41 |
| H1019 | Double-Deck Multi-Function Rack | 25 | 1pc*25bags | 61*41*41 |

**CUVETTES
SAMPLE CUPS**



CUVETTES & SAMPLE CUPS

We imported

optical plastics, and is casted with high precision mould. The inside and outside surface of the product is transparent, level and smooth with superior light transmission performance. It is suitable for all kinds of automatic biochemical analyzer, blood coagulation instrument and spectrophotometer at home and abroad.

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---------------------------------------|----------|------------------------|-----------------|
| A1001 | Tianjin Pulisheng & French stago | 1600 | 50pcs/box*32boxes/ctn | 49*33*33 |
| A1001-1 | Steellex Reaction Vessels | 2304 | 72pcs/box*32boxes/ctn | 49*33*33 |
| A1002-1 | American MD 3000 Biochemical Analysis | 960 | 240pcs/box*4boxes/ctn | 49*33*33 |
| A1005 | BN100 Special Protein Apparatus | 2240 | 70pcs/box*32boxes/ctn | 49*33*33 |
| A1007 | American Pacific, German TECO | 3200 | 100pcs/bag*32bags/ctn | 49*33*33 |
| A1007-1 | MC Coatron 1800 | 3200 | 100pcs/bag*32bags/ctn | 49*33*33 |
| A1008 | Beckman | 1500 | 100pcs/box*15boxes/ctn | 49*33*33 |
| A1008-1 | Beckman Top | 1500 | 100pcs/box*15boxes/ctn | 49*33*33 |

**CUVETTES
SAMPLE CUPS**



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|------------------------------|----------|-------------------|-----------------|
| A1018 | Spain Cruor Apparatus | 8000 | 1000pcs*8bags | 49*33*33 |
| A1019 | German LABOR | 10000 | 1000pcs*10bags | 49*33*33 |
| A1020 | German TECO Sample Cup | 12000 | 1000pcs*12bags | 49*33*33 |
| B1009 | Japan TOA CA1500 | 20000 | 2000pcs*10bags | 48*32*37 |
| B1010 | Japan TOA CA530 | 20000 | 2000pcs*10bags | 48*32*37 |
| B1011 | Japan TOA CA6000 | 20000 | 2000pcs*10bags | 48*32*37 |
| B1020 | Japan TOA CA7000 | 20000 | 1000pcs*20boxes | 48*32*37 |
| B1005 | Abbott AXSYM | 3600 | 400pcs*9bags | 48*33*33 |
| B1016 | For Special Protein Analyzer | 960 | 40pcs/box*24boxes | 49*34*52 |

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|-----------------|-----------------|
| B1025 | Abbott R2000 | 10000 | 1000pcs*10boxes | 48*32*52 |
| B1026 | Siemens ADVIA Centaur CP | 6000 | 1000pcs*6bags | 46*32*29 |
| B1027 | Siphon quantitative extraction tubes | 10000 | 1000pcs*10bags | 48*32*42 |
| B1015 | Rayto RAC-1800 hematology machine | 20000 | 1000pcs*20boxes | 48*32*37 |
| C1001 | Hitachi 7150, 7060 Series; Aulinbas 2ml | 5000 | 500pcs*10bags | 49*33*33 |
| C1002 | BECKMAN COULTER CX Series, 0.5ml | 8000 | 1000pcs*8bags | 49*33*33 |
| C1003 | IL Sample Cup, Monarch Biochemical Analysis Apparatus(0.5ml) | 13000 | 1000pcs*13bags | 49*33*33 |




**CUVETTES
SAMPLE CUPS**



| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|---|----------|----------------|-----------------|
| C1004 | BECKMAN COULTER 700(1.5ml), AmericanDubang, Holland VITALABS SAMPLE CUP | 8000 | 1000pcs*8bags | 49*33*33 |
| C1025 | FC serum sample cup(L) | 10000 | 1000pcs*10bags | 49*33*28 |
| C1026 | FC serum sample cup(S) | 10000 | 1000pcs*10bags | 39*28*27 |
| C1017 | Hitachi Minimum | 6000 | 500pcs*12bags | 49*33*33 |
| C1019 | Hitachi Conical | 9000 | 1000pcs*9bags | 49*33*33 |
| C1020 | Kodak Sample Cup | 12000 | 1000pcs*12bags | 49*33*33 |
| C1028 | Beckman DX1800RV | 12000 | 1000pcs*12bags | 49*33*34 |
| C1029 | Analysis tube | 2400 | 200pcs*12bags | 49*33*37 |




| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|-------------------------------|----------|--------------------|-----------------|
| C1021 | BOWLIN MAN Sample Tube | 40000 | 5000pcs/box*8boxes | 49*33*33 |
| C1021-1 | BOWLIN MAN black pipette tips | 20000 | 1000pcs*20boxes | 49*33*33 |
| C1022 | Kodak Tip | 16000 | 1000pcs*16bags | 49*33*33 |
| C1021-2 | Roche Tip&Cup | 48boxes | 6boxes/bag*8bags | 63*36*23 |
| C1021-5 | Roche Tip (with box) | 120boxes | 120boxes | 43*37*59 |
| C1021-6 | Roche Cup (with box) | 240boxes | 240boxes | 43*39*41 |
| C1021-7 | Roche 801 | 36 boxes | 36 boxes | |



OTHER DISPOSABLE PLASTICWARES

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|-----------------------------------|-----------------|
| H1010 | PP Container with separating plate | 10 | pile up | |
| E1026 | Pipette shelf | 10 | in bulk | 49*33*33 |
| F1019 | Reservoir type "V" | 200 | Individual Packing | 57*35*30 |
| F1023 | 96 deep-well plate, v-bottom | 50 | Individual Packing 50pcs/case | 50*33*26 |
| F1025 | 96 tip comb for deep-well | 100 | Individual Packing 100pcs/case | 52.5*33*28 |
| Q133 | 8OZ Specimen Beaker, W/Temperature label | 500 | 20pcs*25bags | 42*34*34 |
| GDF01 | Disposable face shield,PET,Size:33-22cm | 200 | 1pc/Polyethylene film bags*10pcs/ | |
| | Thickness: 0.18±0.05mm | | PE bags*20PE bags | 67*38.5*39 |
| GDG01 | Toughened Lens to prevent for Scratches, Anti-flog | 100 | 1pc/PEbags*25pcs/box*4boxes | 32*32*30 |
| GDG02 | Eye Glass, Non-sterile | 200 | 10pcs/cs*20cs | 79*34*44 |

OTHER DISPOSABLE PLASTICWARES



PHARMACEUTICAL PACKAGING MATERIALS

| Item No | Description of Goods | Quantity | Packing | Measurement(cm) |
|---------|--|----------|---------------|-----------------|
| AB001 | Specimen Bag 6*9" | 2000 | 100pcs*20bags | 36*36*33 |
| AB002 | Specimen Bag 6*9" | 2000 | 100pcs*20bags | 36*36*33 |
| GD003 | Specimen Bag | | | |
| CP004 | 4ml HDPE Lab Quality Bottle with cap, narrow mouth | 2500 | in bulk | 48*32*26 |
| CP008 | 8ml HDPE Lab Quality Bottle with cap, narrow mouth | 2000 | in bulk | 47*47*42 |
| CP125 | 125ml HDPE Lab Quality Bottle with cap, narrow mouth | 300 | in bulk | 61*41*41 |
| CP250 | 250ml HDPE Lab Quality Bottle with cap, narrow mouth | 200 | in bulk | 63*44*45 |

